

Anti-SLPI Antibody (9W363)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	9W363
Purification:	Protein A

Applications

Application:	ELISA
Recommended	ELISA: 1:5000-1:10000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human SLPI Protein (TMPY-02389)
Antigen Species:	Human
Synonyms:	WFDC4;ALP;HUSI;WAP4;SLPI;MPI;HUSI-1;BLPI;ALK1
Biology Area:	Serine Proteases and Regulators

Research Background

Secretory leukoprotease inhibitor (SLPI), also called antileukoprotease (ALP), is a 12-kDa, nonglycosylated serine protease inhibitor present in mucous secretions. It is thought to play a role in protecting the mucosae from injury associated with inflammation. SLPI is locally produced by serous cells, including bronchial submucosal glands. Elafin and SLPI are members of larger families of proteins secreted predominantly at mucosal sites, and have been shown to be modulated in multiple pathological conditions. Elafin and SLPI are structurally related in that both have a fold with a four-disulfide core or whey acidic protein (WAP) domain responsible for inhibiting proteases. SLPI is a prominent innate immune protein of the respiratory tract, possessing serine protease inhibitor activity, antibacterial activity, and anti-inflammatory/immunomodulatory activity.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481