

## Anti-Hsp27/HSPB1 Antibody-HRP (3Y390)

### Product Details

Ig Type:	Mouse IgG2b
Reactivity:	Human
Conjugation:	HRP
Clone:	3Y390
Purification:	Protein A

### Applications

Application:	ELISA
Recommended	ELISA: 0.1-1 µg/ml

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free. Keep away from direct sunlight.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Human HSP27 Protein (TMPY-00823)
Antigen Species:	Human
Synonyms:	HMN2B;HSP28;Hsp27;CMT2F;HSPB1;SRP27;HSP 27
Biology Area:	Cancer Drug Targets

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481