

Anti-Hsp27/HSPB1 Antibody-PE (20259)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	PE
Clone:	20259
Purification:	Protein A

Applications

Verified Activity:	Flow cytometric analysis of Human HSP27 expression in HeLa cells. The cells were treated according to manufacturer's manual (BD Pharmingen™ Cat. No. 554714), and then stained with PE Rabbit anti-Human HSP27. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of intact cells.
Application:	FCM
Recommended	5 µl/Test, 0.1 mg/ml

Properties

Stability & Storage:	Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight. Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human HSP27 / HSPB1 protein (TMPY-00823)
Antigen Species:	Human
Synonyms:	HMN2B;SRP27;HSP 27;HSP28;HSPB1;Hsp27;CMT2F
Biology Area:	Cancer Drug Targets

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481