

Anti-IL-33 Antibody-FITC (4M737)

Product Details

Ig Type:	Mouse IgG2b
Reactivity:	Human
Conjugation:	FITC
Clone:	4M737
Purification:	Protein A

Applications

Verified Activity:	Flow cytometric analysis of Human IL-33 on HUVEC cells. HUVEC cells were treated according to manufacturer's manual (BD Pharmingen™ Cat. No. 554714), and stained with FITC Mouse anti-IL-33 (Bold line hisgram).
Application:	FCM
Recommended	10 µl/Test, 0.1 mg/ml

Properties

Stability & Storage:	Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight. Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human IL-33 protein (TMPY-00861)
Antigen Species:	Human

Research Background

Interleukin 33 (IL-33), also known as DVS27 or NF-HEV (Nuclear Factor from High Endothelial Venules), is a pro-inflammatory protein and a chromatin-associated cytokine of the IL-1 family with high sequence and structural similarity to IL-1 and IL-18. IL-33 protein is expressed highly and rather selectively by high endothelial venule endothelial cells (HEVECs) in human tonsils, Peyer's patches, and lymph nodes. IL-33 protein has transcriptional regulatory properties, and the researches suggested that IL-33 is a dual-function protein that might act both as a cytokine and as an intracellular nuclear factor. As a type 2 cytokines, IL-33 protein also play a pivotal role in helminthic infection and allergic disorders.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481