

## Anti-CD32/FCGR2 Antibody (3K55)

### Product Details

Ig Type:	Mouse IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	3K55
Purification:	Protein A

### Applications

Verified Activity:	Flow cytometric analysis of Human CD32 expression on human whole blood lymphocytes (left panel), monocytes (middle panel) and granulocytes (right panel). Cells were stained with purified anti-Human CD32, then a FITC-conjugated second step antibody. The histograms were derived from gated events with the forward and side light-scatter characteristics of viable lymphocytes, monocytes and granulocytes.
Application:	ELISA,FCM
Recommended	ELISA: 1:1000-1:2000; FCM: 1:25-1:100

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Human CD32a protein (TMPY-00764)
Antigen Species:	Human
Synonyms:	FCG2;FcgammaRIIa;CDw32;CD32A;FCGR2;CD32;IGFR2;FcGR;Fc-gamma-RIIa;FCGR2A1
Biology Area:	Fc Receptors, ITIM/ITAM Immunoreceptors and Related Molecules

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481