

Anti-CD32/FCGR2 Antibody-PE (3K55)

Product Details

Ig Type:	Mouse IgG1
Reactivity:	Human
Conjugation:	PE
Clone:	3K55
Purification:	Protein A

Applications

Verified Activity:	Flow cytometric analysis of Human CD32 expression on human whole blood lymphocytes (left panel), monocytes (middle panel) and granulocytes (right panel). Cells were stained with PE-conjugated anti-Human CD32. The histograms were derived from gated events with the forward and side light-scatter characteristics of viable lymphocytes, monocytes and granulocytes.
Application:	FCM
Recommended	10 µl/Test, 0.1 mg/ml

Properties

Stability & Storage:	Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight. Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human CD32a/FCGR2A/Fc gamma RIIA Protein (TMPY-00764)
Antigen Species:	Human
Synonyms:	Fc-gamma-RIIa; CD32; CDw32; FCGR2; FCGR2A1; FCG2; CD32A; Fc gamma RIIa; FcGR; IGFR2
Biology Area:	Fc Receptors, ITIM/ITAM Immunoreceptors and Related Molecules

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481