

## Anti-CD32A Antibody-FITC (2A526)

### Product Details

Ig Type:	Mouse IgG1
Reactivity:	Human
Conjugation:	FITC
Clone:	2A526
Purification:	Protein A

### Applications

Verified Activity:	Flow cytometric analysis of Human CD32a expression on human whole blood lymphocytes (left panel), monocytes (middle panel) and granulocytes (right panel). Cells were stained with FITC-conjugated anti-Human CD32a. The histograms were derived from gated events with the forward and side light-scatter characteristics of viable lymphocytes, monocytes and granulocytes.
Application:	FCM
Recommended	10 µl/Test, 0.1 mg/ml

### Properties

Stability & Storage:	Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Human CD32a/FCGR2A/Fc gamma RIIA Protein (TMPY-00764)
Antigen Species:	Human
Synonyms:	Fc fragment of IgG receptor IIa
Biology Area:	Fc Receptors, ITIM/ITAM Immunoreceptors and Related Molecules

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481