

Anti-MMP-3 Antibody (3U852)

Product Details

Ig Type:	Mouse IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	3U852
Purification:	Protein A

Applications

Application:	ELISA,ELISA(Det)
Recommended	ELISA: 1:1000-1:2000; ELISA(Det): 1:1000-1:10000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human MMP-3 protein (TMPY-02689)
Antigen Species:	Human
Synonyms:	STMV;SL-1;CHDS6;STR1;matrix metalloproteinase 3;MMP-3;STMV1

Research Background

Matrix metalloproteinase 3 (abbreviated as MMP3) is also known as stromelysin 1 and progelatinase. MMP3 is a member of the matrix metalloproteinase (MMP) family whose members are involved in the breakdown of extracellular matrix in normal physiological processes, such as embryonic development, reproduction, tissue remodeling, and disease processes including arthritis and metastasis. As a secreted zinc-dependent endopeptidase, MMP3 exerts its functions mainly in the extracellular matrix. This protein is activated by two major endogenous inhibitors: alpha2-macroglobulin and tissue inhibitors of metalloproteinases (TIMPs). MMP3 plays a central role in degrading collagen types II, III, IV, IX, and X, proteoglycans, fibronectin, laminin, and elastin. Also, MMP3 can activate other MMPs such as MMP1, MMP7, and MMP9, rendering MMP3 crucial in connective tissue remodeling. Dysregulation of MMPs has been implicated in many diseases including arthritis, chronic ulcers, encephalomyelitis, and cancer. Synthetic or natural inhibitors of MMPs result in inhibition of metastasis, while up-regulation of MMPs led to enhanced cancer cell invasion.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481