

Anti-His Tag Antibody (3U153)

Product Details

Ig Type:	Mouse IgG2b
Reactivity:	other
Conjugation:	Unconjugated
Clone:	3U153
Purification:	Protein A

Applications

Application:	ELISA
Recommended	ELISA: 1:5000-1:20000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthetic peptide: 6×His
Antigen Species:	Other

Research Background

The hexahistidine tag (6× His-tag) is the most frequently used affinity tag for protein enrichment. His-tagged proteins can be purified easily by the chelated metal ions as affinity ligands. The basis for affinity purification is known as immobilized metal affinity chromatography (IMAC). His-tag can bind best to IMAC resin in near-neutral buffer conditions (physiologic pH and ionic strength), and thus the fusion proteins can be eluted with binding buffer containing certain concentrations of imidazole. If possible, the elution is also accomplished with low pH (e.g., 0.1M glycine-HCl, pH 2.5) or an excess of strong chelator (e.g., EDTA).

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481