

Anti-TFPI Antibody (6I399)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	6I399
Purification:	Protein A

Applications

Verified Activity:	<p>1. Anti-TFPI rabbit monoclonal antibody at 1:500 dilution.</p> <ul style="list-style-type: none">-Lane A: MCF7 Whole Cell lysate.-Lysates/proteins at 30 µg per lane.-Secondary-Goat Anti-Rabbit IgG H&L (Dylight800) at 1/10000 dilution.-Developed using the Odyssey technique.-Performed under reducing conditions.-Predicted band size:35 kDa.-Observed band size:40 kDa. <p>2. TFPI was immunoprecipitated using:</p> <ul style="list-style-type: none">-Lane A:0.5 mg MCF-7 Whole Cell Lysate-0.5 µL anti-TFPI rabbit monoclonal antibody and 15 µl of 50 % Protein G agarose.-Primary antibody:-Anti-TFPI rabbit monoclonal antibody, at 1:1000 dilution.-Secondary antibody:-Dylight 800-labeled antibody to rabbit IgG (H+L), at 1:5000 dilution.-Developed using the odyssey technique.-Performed under reducing conditions.-Predicted band size: 35 kDa.-Observed band size: 45 kDa
Application:	IP,WB
Recommended	WB: 1:500-1:1000; IP: 0.2-1 µL/mg of lysate

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human TFPI Protein (TMPY-01030)
Antigen Species:	Human
Synonyms:	EPI;LACI;TFI;TFPI1;tissue factor pathway inhibitor (lipoprotein-associated coagulation inhibitor)
Biology Area:	Serine Proteases and Regulators

Research Background

Tissue factor pathway inhibitor (TFPI) is the natural inhibitor of TF coagulant and signaling activities. It is a Kunitz-type serine proteinase inhibitor that down-regulates tissue factor-initiated blood coagulation. With its Kunitz domains, TFPI exhibits significant homology with human inter-alpha-trypsin inhibitor and bovin basic pancreatic trypsin inhibitor. TFPI is the natural inhibitor of TF coagulant and signaling activities. The importance of TFPI in the regulation of blood coagulation is emphasized by how its activity is modulated in human disease. In a factor (F) Xa-dependent feedback system, TFPI binds directly and inhibits the TF-FVII/FVIIa complex. Normally, TFPI exists in plasma both as a full-length molecule and as variably carboxy-terminal truncated forms. TFPI also circulates in complex with plasma lipoproteins. The levels and the dual inhibitor effect of TFPI on FXa and TF-FVII/FVIIa complex offers insight into the mechanisms of various pathological conditions triggered by TF. TFPI may play an important role in modulating TF-induced thrombogenesis and it may also provide a unique therapeutic approach for prophylaxis and/or treatment of various diseases. In addition, studies have shown that TFPI exhibits antiangiogenic and antimetastatic effects in vitro and in vivo. In animal models of experimental metastasis, both circulating and tumor cell-associated TFPI are shown to significantly reduce tumor cell-induced coagulation activation and lung metastasis.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481