

Anti-BMP-4 Antibody (8A124)

Product Details

Ig Type:	Mouse IgG2a
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	8A124
Purification:	Protein A

Applications

Application:	ELISA
Recommended	ELISA: 1:1000-1:2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human BMP4 Protein (TMPY-06842)
Antigen Species:	Human
Synonyms:	MCOPS6;OFC11;ZYME;BMP2B;BMP2B1;bone morphogenetic protein 4
Biology Area:	Early Mesodermal Lineage Markers, Cardiac Stem Cell Markers

Research Background

BMP4 (bone morphogenetic protein-4)-directed differentiation of human embryonic stem cells (ESCd). Autophagy is a conserved catabolic process with complicated roles in tumor development. Bone morphogenetic protein 4 (BMP4), a member of the transforming growth factor (TGF- β) family of regulatory proteins, plays a crucial role in human malignancies. BMP4 treatment promoted HCC cells proliferation and induced autophagy. Mechanistic study revealed that the induction of autophagy by BMP4 was mediated through activating the JNK1/Bcl2 pathway. And the JNK1 inhibitor and knockdown of JNK1 could attenuate autophagy induced by BMP4 and eliminated BMP4-promoted HCC cells growth. BMP4 promoted HCC proliferation by autophagy activation through JNK1/Bcl-2 signaling.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481