

Anti-CCL28 Antibody (6R663)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	6R663
Purification:	Protein A

Applications

Application:	ELISA
Recommended	ELISA: 1:5000-1:10000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human CCL28 / SCYA28 Protein (TMPY-03728)
Antigen Species:	Human
Synonyms:	MEC;SCYA28;CCK1;chemokine (C-C motif) ligand 28

Research Background

CCL28 chemokine is expressed by epithelial cells of various mucosal tissues. This chemokine binds to CCR3 and CCR10 receptors and plays an essential role in the IgA antibody secreting cells (IgA-ASC) homing to mucosal surfaces and lactating mammary gland as well. Besides, CCL28 has been shown to exert a potent antimicrobial activity against both Gram-negative and Gram-positive bacteria and fungi. The potent antimicrobial function of CCL28 combined with its wide distribution in mucosal tissues and secretions suggest that this protein plays an important role in innate immune protection of the epithelial surfaces. CCL28 is a human chemokine constitutively expressed by epithelial cells in diverse mucosal tissues and is known to attract a variety of immune cell types including T-cell subsets and eosinophils. Elevated levels of CCL28 have been found in the airways of individuals with asthma, and previous studies have indicated that CCL28 plays a vital role in the acute development of post-viral asthma. CCL28 presents a novel target for the development of alternative asthma therapeutics. The dental decay of children leads to the secretion of chemokine CCL28, which promotes the secretion of sIgA in saliva. CC chemokine ligand28 (CCL28) has been reported as a severity marker in atopic dermatitis.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481