

Anti-MPIF-1/CCL23 Antibody (5S851)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	5S851
Purification:	Protein A

Applications

Application:	ELISA(Cap)
--------------	------------

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human CCL23 / MIP 3 Protein
Antigen Species:	Human
Synonyms:	Ckb-8-1;CK-BETA-8;chemokine (C-C motif) ligand 23;MPIF-1;MPIF1;CCL23;MIP-3;SCYA23;Ckb-8;MIP3;CK-β-8;CKb8;hmrp-2a

Research Background

CCL23, also known as MIP 3, is a small cytokine which belongs to the CC chemokine family. Cytokines are a family of secreted proteins involved in immunoregulatory and inflammatory processes. CCL23 is predominantly expressed in lung and liver tissue, but also can be detected in bone marrow and placenta. It displays chemotactic activity on resting T lymphocytes and monocytes, lower activity on neutrophils and no activity on activated T lymphocytes. CCL23 is also a strong suppressor of colony formation by a multipotential hematopoietic progenitor cell line.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481