

## Anti-MICB Antibody (2D186)

## Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	2D186
Purification:	Protein A

## Applications

	Anti-MICB rabbit monoclonal antibody at 1:500 dilution. -Lane A: U87MG Whole Cell lysate. -Lysates/proteins at 30 µg per lane. -Secondary
Verified Activity:	-Goat Anti-Rabbit IgG H&L (Dylight800) at 1/10000 dilution. -Developed using the Odyssey technique. -Performed under reducing conditions. -Predicted band size:40 kDa. -Observed band size:67 kDa
Application:	WB
Recommended	WB: 1:500-1:1000

## Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

## Antigen Details

Immunogen:	Recombinant Protein: Human MICB protein (TMPY-01066)
Antigen Species:	Human
Synonyms:	MHC class I polypeptide-related sequence B;PERB11.2

## Research Background

MHC class I polypeptide-related sequence B, also known as MICB, is a heavily glycosylated protein serving as a ligand for the type II receptor NKG2D. MICB shares 85% amino acid identity with MICA, a closely related protein, both of which contain three extracellular immunoglobulin-like domains, but without the capacity to bind peptide or interact with beta-2-microglobulin. acting as a stress-induced self-antigen, binding of MICB to the NKG2D receptor activates the cytolytic response of natural killer (NK) cells, CD8+αβ T cells, and γδ T cells on which the receptor is expressed. MICA/B is minimally expressed on normal cells, but are frequently expressed on epithelial tumors and can be induced by bacterial and viral infections. MICA/B recognition thus is involved in tumor surveillance, viral infections, and autoimmune diseases.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481