

Anti-CD40 Antibody-APC (8Z859)

Product Details

Ig Type:	Mouse IgG2b
Reactivity:	Human
Conjugation:	APC
Clone:	8Z859
Purification:	Protein A

Applications

Verified Activity:	Flow cytometric analysis of human CD40 expression on human whole blood lymphocytes. Human whole blood lymphocytes were stained with APC-conjugated anti-Human CD40 and PerCP-Cy™5.5 conjugated anti-Human CD19 (BD Pharmingen™ Cat. No. 561295). The dot plots were derived from events with the forward and side light-scatter characteristics of viable lymphocytes.
Application:	FCM
Recommended	5 µl/Test, 0.1 mg/ml

Properties

Stability & Storage:	Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight. Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human CD40 / TNFRSF5 protein (TMPY-01181)
Antigen Species:	Human
Synonyms:	TRAP;CD40 molecule, TNF receptor superfamily member 5;T-BAM;p50;Tnfrsf5;AI326936;Igm;gp39;IMD3;HIGM1;Bp50
Biology Area:	Cancer Drug Targets

Research Background

CD40, also known as TNFRSF5, is a member of the TNF receptor superfamily which are single transmembrane-spanning glycoproteins. CD40 protein plays an essential role in mediating a broad variety of immune and inflammatory responses including T cell-dependent immunoglobulin class switching, memory B cell development, and germinal center formation. CD40 protein is expressed in B cells, dendritic cells, macrophages, endothelial cells, and several tumor cell lines. Defects in CD40 result in hyper-IgM immunodeficiency type 3 (HIGM3). In addition, CD40/CD40L interaction is found to be necessary for amyloid-beta-induced microglial activation, and thus is thought to be an early event in Alzheimer disease pathogenesis. Cancer Immunotherapy Co-stimulatory Immune Checkpoint Targets Immune Checkpoint Immune Checkpoint Detection: Antibodies Immune Checkpoint Detection: ELISA Antibodies Immune Checkpoint Detection: FCM Antibodies Immune Checkpoint Detection: ICC Antibodies Immune Checkpoint Detection: IP Antibodies Immune Checkpoint Detection: WB Antibodies Immune Checkpoint Proteins Immune Checkpoint Targets Immunotherapy Targeted Therapy

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481