

Anti-14-3-3 epsilon Antibody (2S60)

Product Details

Ig Type:	Mouse IgG2a
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	2S60
Purification:	Protein A

Applications

Verified Activity:	<p>1. Anti-YWHAE mouse monoclonal antibody at 1:500 dilution.</p> <ul style="list-style-type: none">-Lane A: HepG2 Whole Cell Lysate.-Lane B: Jurkat Whole Cell Lysate.-Lane C: U-251 MG Whole Cell Lysate.-Lane D: A375 Whole Cell lysate. <p>-Lysates/proteins at 30 µg per lane.</p> <p>-Secondary</p> <ul style="list-style-type: none">-Goat Anti-Mouse IgG (H+L)/HRP at 1/10000 dilution.-Developed using the ECL technique.-Performed under reducing conditions.-Predicted band size:29 kDa.-Observed band size:33 kDa. <p>2. YWHAE was immunoprecipitated using:</p> <ul style="list-style-type: none">-Lane A:0.5 mg Jurkat Whole Cell Lysate.-Lane B:0.5 mg HepG2 Whole Cell Lysate.-Lane C:0.5 mg SH-SY5Y Whole Cell Lysate. <p>-4 µL anti-YWHAE mouse monoclonal antibody and 60 µg of Immunomagnetic beads Protein A/G.</p> <p>-Primary antibody:</p> <ul style="list-style-type: none">-Anti-YWHAE mouse monoclonal antibody, at 1:100 dilution. <p>-Secondary antibody:</p> <ul style="list-style-type: none">-Rabbit Anti-Mouse IgG F(ab)₂/HRP at 1/10000 dilution.-Developed using the ECL technique.-Performed under reducing conditions.-Predicted band size: 34 kDa.-Observed band size:34 kDa
Application:	IP,WB

A DRUG SCREENING EXPERT

Recommended WB: 1:500-1:2000; IP: 1-5 µL/mg of lysate

Properties

Stability & Storage: Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.

Shipping: Shipping with blue ice.

Antigen Details

Immunogen: Recombinant Protein: Human 14-3-3 epsilon Protein (TMPY-01274)

Antigen Species: Human

Synonyms: 14-3-3 ε; tyrosine 3-monooxygenase/tryptophan 5-monooxygenase activation protein, ε; tyrosine 3-monooxygenase/tryptophan 5-monooxygenase activation protein, epsilon; AU019196

Biology Area: phospho-Serine/phospho-Threonine Binding Proteins, Adaptor Proteins, Apoptosis Adaptor Proteins

Research Background

YWHAE, also known as 14-3-3 epsilon, mediate signal transduction by binding to phosphoserine-containing proteins. 14-3-3 epsilon / YWHAE is a member of the 14-3-3 proteins family. 14-3-3 proteins are a group of highly conserved proteins that are involved in many vital cellular processes such as metabolism, protein trafficking, signal transduction, apoptosis and cell cycle regulation. 14-3-3 proteins are mainly localized in the synapses and neuronal cytoplasm, and seven isoforms have been identified in mammals. This family of proteins was initially identified as adaptor proteins which bind to phosphoserine-containing motifs. Binding motifs and potential functions of 14-3-3 proteins are now recognized to have a wide range of functional relevance. 14-3-3 epsilon / YWHAE is found in both plants and mammals, and this protein is 100% identical to the mouse ortholog. YWHAE interacts with CDC25 phosphatases, RAF1 and IRS1 proteins, suggesting its role in diverse biochemical activities related to signal transduction, such as cell division and regulation of insulin sensitivity. It has also been implicated in the pathogenesis of small cell lung cancer. 14-3-3 epsilon / YWHAE is implicated in the regulation of a large spectrum of both general and specialized signaling pathways. 14-3-3 epsilon / YWHAE Binds to a large number of partners, usually by recognition of a phosphoserine or phosphothreonine motif. This Binding generally results in the modulation of the activity of the binding partner.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481
