

Anti-IL-17RA Antibody (3M821)

Product Details

Ig Type:	Mouse IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	3M821
Purification:	Protein A

Applications

	Anti-IL17RA mouse monoclonal antibody at 1:500 dilution. -Lane A: Hela Whole Cell lysate. -Lysates/proteins at 30 µg per lane. -Secondary
Verified Activity:	-Goat Anti-Mouse IgG H&L (Dylight800) at 1/15000 dilution. -Developed using the Odyssey technique. -Performed under reducing conditions. -Predicted band size:96 kDa. -Observed band size:105 kDa(We are unsure as to the identity of these extra bands.)
Application:	ELISA,ELISA(Cap),WB
Recommended	WB: 1:500-1:1000; ELISA: 1:1000-1:2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human IL17RA protein (TMPY-01371)
Antigen Species:	Human
Synonyms:	hIL-17R;CANDF5;IL-17RA;interleukin 17 receptor A;CDw217;IL17R;CD217

Research Background

Interleukin-17 receptor (IL-17R), also known as Interleukin-17 receptor A (IL-17RA) and CD217 antigen (CD217), is a cytokine receptor that binds interleukin 17. IL-17R/IL-17RA (CD217) is a proinflammatory cytokine secreted by activated T-lymphocytes. It is a potent inducer of the maturation of CD34-positive hematopoietic precursors into neutrophils. IL-17R/IL-17RA (CD217) is a ubiquitous type I membrane glycoprotein that binds with low affinity to interleukin 17A. Interleukin 17A and its receptor IL-17RA play a pathogenic role in many inflammatory and autoimmune diseases such as rheumatoid arthritis. Like other cytokine receptors, this receptor likely has a multimeric structure. Defects in IL-17R/IL-17RA (CD217) are the cause of familial candidiasis type 5 (CANDF5). CANDF5 is a rare disorder with altered immune responses and impaired clearance of fungal infections, selective against *Candida*. It is characterized by persistent and/or recurrent infections of the skin, nails, and mucous membranes caused by organisms of the genus *Candida*, mainly *Candida albicans*.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481