

## Anti-Flag Tag Antibody-HRP (6Z668)

### Product Details

Ig Type:	Mouse IgG2a
Conjugation:	HRP
Clone:	6Z668
Purification:	Protein A

### Applications

Application:	ELISA
Recommended	ELISA: 0.1-1 µg/mL

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free. Keep away from direct sunlight.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthetic peptide: DYKDDDDK-tag sequence.
Antigen Species:	Other

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481