

## Anti-Cadherin 16/CDH16 Antibody (9V223)

### Product Details

Ig Type:	Mouse IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	9V223
Purification:	Protein A

### Applications

Verified Activity:	Immunochemical staining of human CDH16 in human kidney with mouse monoclonal antibody (1:60, formalin-fixed paraffin embedded sections).
Application:	IHC-P
Recommended	IHC-P: 1:50-1:200

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Human KSP-Cadherin / Cadherin-16 Protein (TMPY-01763)
Antigen Species:	Human
Synonyms:	cadherin 16;UNQ695/PRO1340

### Research Background

KSP-Cadherin/Cadherin-16 is a member of the cadherin superfamily, calcium-dependent, membrane-associated glycoproteins. The protein consists of an extracellular domain containing 6 cadherin domains, a transmembrane region and a truncated cytoplasmic domain but lacks the prosequence and tripeptide HAV adhesion recognition sequence typical of most classical cadherins. Expression is exclusively in kidney, where the protein functions as the principal mediator of homotypic cellular recognition, playing a role in the morphogenic direction of tissue development. KSP-Cadherin/Cadherin-16 can be detected at later stages of tubulogenesis during human renal development and in the distal tubules of adult kidneys, no expression was found by immunohistochemistry or Western blot analysis in RCC tumour tissues and several RCC cell lines. However, despite the lack of protein expression, mRNA synthesis of KSP-Cadherin/Cadherin-16 could be detected by reverse transcriptase-polymerase chain reaction analysis in all RCC tissues and most of the RCC cell lines studied, although at a reduced level. The loss of KSP-Cadherin/Cadherin-16 protein was only observed in the malignant part of the tumour kidneys, whereas in the normal part of the affected kidneys KSP-Cadherin/Cadherin-16 expression was clearly detected. These results indicate a downregulation of Ksp-cadherin in RCC and suggest a role for this cell adhesion molecule in tumour suppression.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481