

Anti-Neuronostatin Antibody (2M309)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	2M309
Purification:	Protein A

Applications

Application:	ELISA,ELISA(Cap)
Recommended	ELISA: 1:5000-1:10000; ELISA(Cap): 1:250-1:2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthetic peptide: center region of Human Neuronostatin
Antigen Species:	Human
Synonyms:	SST1;Smst

Research Background

Somatostatin (SST) is a cyclic peptide that is understood to inhibit the release of hormones and neurotransmitters from a variety of cells by binding to one of five canonical G protein-coupled SST receptors (SSTR1 to SSTR5). Somatostatin (SST) is an inhibitor of the release of somatotropin, glucagon, insulin, gastrointestinal hormones, and other secretory proteins (1). SST is also known as somatotropin release-inhibiting factor (SRIF). SST is a cyclic polypeptide with two biologically active isoforms, SRIF-14 and SRIF-28, of 14 and 28 amino acids, respectively. SST negatively regulates cell proliferation and the release of multiple hormones via activation of its cognate receptors. A variety of SST analogs, some with high affinity and selectivity of receptor subtypes, have been synthesized and developed.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481