

## Anti-GLP-1/Glucagon-like peptide 1 Antibody (8P428)

### Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	8P428
Purification:	Protein A

### Applications

Application:	ELISA,ELISA(Cap)
Recommended	ELISA: 1:5000-1:10000; ELISA(Cap): 1:250-1:2000

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthetic peptide: Human GLP-1/Glucagon-like
Antigen Species:	Human
Synonyms:	GRPP;GLP1

### Research Background

Variants in GCG may associate with type 2 diabetes, obesity and/or related metabolic traits. WNT signaling pathway in a pivotal role in regulating pancreas development as well as islet function, insulin production and secretion, GCG is one of the key effectors in the WNT signaling pathway.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481