

Anti-N Cadherin Antibody (3E575)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	3E575
Purification:	Protein A

Applications

1. Immunofluorescence staining of Human CDH2 in HeLa cells. Cells were fixed with 4% PFA, blocked with 10% serum, and incubated with rabbit anti-Human CDH2 monoclonal antibody (1:60). Then cells were stained with the Alexa Fluor® 488-conjugated Goat Anti-rabbit IgG secondary antibody (green) and counterstained with DAPI (blue). Positive staining was localized to membrane.

2. Flow cytometry analysis of N-Cadherin on HeLa cells. HeLa cells were harvested with (Left) or without (Right) trypsinization [please note, the epitope is sensitive to trypsin] and stained with either purified Rabbit IgG isotype control (dashed line) or Rabbit Anti-Human CD325 antibody [solid line, and Clone 1T666 (upper), Clone 3E575 (bottom)] at matching concentrations, then a FITC-conjugated second step antibody.

3. Anti-CDH2 rabbit monoclonal antibody at 1:500 dilution.

- Lane A: HeLa Whole Cell lysate.
- Lysates/proteins at 30 µg per lane.
- Secondary

Verified Activity: -Goat Anti-Rabbit IgG H&L (Dylight800) at 1/10000 dilution.

- Developed using the Odyssey technique.
- Performed under reducing conditions.
- Predicted band size:99 kDa.
- Observed band size:130 kDa.

4. Anti-N Cadherin rabbit monoclonal antibody at 1:500 dilution.

- Lane A: 293T Whole Cell Lysate.
- Lane B: N Cadherin konckout 293T Whole Cell lysate.
- Lysates/proteins at 30 µg per lane.
- Secondary
- Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution.
- Developed using the ECL technique.
- Performed under reducing conditions.
- Predicted band size:99 kDa.
- Observed band size:130 kDa(Validation Experiment)

A DRUG SCREENING EXPERT

Application: ELISA,FCM,ICC/IF,WB

Recommended WB: 1:500-1:2000; ELISA: 1:5000-1:10000; ICC-IF: 1:20-1:100; FCM: 1:25-1:100

Properties

Stability & Storage: Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.

Shipping: Shipping with blue ice.

Antigen Details

Immunogen: Recombinant Protein: Human N-Cadherin / CD325 / CDH2 Protein (TMPY-01143)

Antigen Species: Human

Synonyms: N-cadherin;CDHN;cadherin 2, type 1, N-cadherin (neuronal);Ncad

Biology Area: Early Mesodermal Lineage Markers

Research Background

Cadherins are calcium-dependent cell adhesion proteins, and they preferentially interact with themselves in a homophilic manner in connecting cells. Cadherin 2 (CDH2), also known as N-Cadherin (neuronal) (NCAD), is a single-pass transmembrane protein and a cadherin containing 5 cadherin domains. N-Cadherin displays a ubiquitous expression pattern but with different expression levels between endocrine cell types. CDH2 (NCAD) has been shown to play an essential role in normal neuronal development, which is implicated in an array of processes including neuronal differentiation and migration, and axon growth and fasciculation. In addition, N-Cadherin expression was upregulated in human HSC during activation in culture, and function or expression blocking of N-Cadherin promoted apoptosis. During apoptosis, N-Cadherin was cleaved into 20-100 kDa fragments. It may provide a novel target for therapies that are directed toward intimal proliferative disorders, including restenosis and vascular bypass graft failure. N-Cadherin is associated with tumor aggressiveness and metastatic potential and may contribute to tumor progression.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481
