

Anti-CD16b Antibody (7T676)

Product Details

Ig Type:	Mouse IgG2a
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	7T676
Purification:	Protein A

Applications

Application:	ELISA
Recommended	ELISA: 1:1000-1:2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human CD16b protein (TMPY-00701)
Antigen Species:	Human
Synonyms:	CD16b;FCR-10;Fc gamma RIIIb;FCRIIIb;Fc fragment of IgG receptor IIIb;CD16;Fc γ RIIIb;FCRIII;FCGR3;FCG3
Biology Area:	Fc Receptors, ITIM/ITAM Immunoreceptors and Related Molecules

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481