

## Anti-CD16b Antibody (3W146)

### Product Details

Ig Type:	Mouse IgG2a
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	3W146
Purification:	Protein A

### Applications

Verified Activity:	1. Immunochemical staining of human CD16b in human liver with mouse monoclonal antibody at 1:200 dilution, formalin-fixed paraffin embedded sections.
	2. Immunochemical staining of human CD16b in human spleen with mouse monoclonal antibody at 1:200 dilution, formalin-fixed paraffin embedded sections.
	3. Immunochemical staining of human CD16b in human lung with mouse monoclonal antibody at 1:200 dilution, formalin-fixed paraffin embedded sections.
Application:	IHC-P
Recommended	IHC-P: 1:100-1:1000

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Human CD16b / FCGR3B protein (TMPY-00701)
Antigen Species:	Human
Synonyms:	FCRIII;Fc gamma RIIB;CD16b;Fc fragment of IgG receptor IIIB;FCGR3;Fc gamma RIIB;FCR-10;CD16;FCRIIB;FCG3
Biology Area:	Fc Receptors, ITIM/ITAM Immunoreceptors and Related Molecules

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481