

## Anti-Phospho-MAPKAPK-2 (Thr334) Antibody (5R815)

## Product Details

Ig Type:	Rabbit IgG Human;
Reactivity:	Predicted to React with:Species predicted to react based on 100% sequence homology: Mouse, Rat, Bovine, Rabbit
Conjugation:	Unconjugated
Clone:	5R815
Purification:	Protein A

## Applications

Verified Activity:	<ol style="list-style-type: none"><li>Western blot analysis of extracts from serum-starved A431, untreated(line A) or treated with calyculin A (100mM, 30min; +)(line B), using Phospho-MAPKAPK-2 (Thr334) rabbit monoclonal Antibody at 1:200 dilution (upper) or Anti-Actin Antibody, Chimeric Rabbit Monoclonal at 1:50000 dilution (lower).</li><li>Western blot analysis of extracts from serum-starved A431, untreated (line A); treated with calyculin A (100mM, 30min; +) (line B); treated with calyculin A and <math>\lambda</math>-phosphatase (line C) using Phospho-MAPKAPK-2 (Thr334) rabbit monoclonal Antibody at 1:200 dilution. (Validation Experiment)</li><li>Western blot analysis of extracts from serum-starved A431, untreated (line A); treated with calyculin A (100mM, 30min), without peptide (line B) or antigen-specific phosphopeptide (line C) or antigen-specific peptide (line D) using Phospho-MAPKAPK-2 (Thr334) rabbit monoclonal Antibody at 1:200 dilution. (Validation Experiment)</li></ol>
Application:	WB
Recommended	WB: 1:500-1:2000

## Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

## Antigen Details

Immunogen:	A synthetic peptide: residues around (Thr334) of Human Phospho-MAPKAPK-2
Antigen Species:	Human
Synonyms:	Phospho-MAPKAPK-2 (T334);MAPKAPK-2;MK2;MK-2;MAPKAPK-2 (p-T334);MAPKAPK-2 (p-Thr334);MAPKAPK2;MAPKAP-K2;p-MAPKAPK-2 (T334);p-MAPKAPK-2 (Thr334)

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481