

Anti-CD68 Antibody-PE (9Z25)

Product Details

| | |
|---------------|------------|
| Ig Type: | Mouse IgG1 |
| Reactivity: | Human |
| Conjugation: | PE |
| Clone: | 9Z25 |
| Purification: | Protein A |

Applications

| | |
|--------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Verified Activity: | Flow cytometric analysis of human CD68 expression on human whole blood monocytes. The cells were treated according to manufacturer's manual (BD Pharmingen™ Cat. No. 554714), stained with PE-conjugated anti-Human CD68. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of viable monocytes. |
| Application: | FCM |
| Recommended | 10 µl/Test, 0.1 mg/ml |

Properties

| | |
|----------------------|--------------------------------------------------------------------------------|
| Stability & Storage: | Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight. |
| Shipping: | Shipping with blue ice. |

Antigen Details

| | |
|------------------|------------------------------------------------------|
| Immunogen: | Recombinant Protein: Human CD68 Protein (TMPY-02972) |
| Antigen Species: | Human |
| Synonyms: | LAMP4;GP110;SCARD1;CD68 molecule |

Research Background

Macrosialin, also known as CD68 and Gp11, is a single-pass type I membrane protein which belongs to the LAMP family. CD68 is highly expressed by blood monocytes and tissue macrophages. It is also expressed in lymphocytes, fibroblasts and endothelial cells. CD68 is expressed in many tumor cell lines which could allow them to attach to selectins on vascular endothelium, facilitating their dissemination to secondary sites. CD68 plays a role in phagocytic activities of tissue macrophages, both in intracellular lysosomal metabolism and extracellular cell-cell and cell-pathogen interactions. It is a commonly used marker for macrophages. However, a number of studies have shown that CD68 antibodies react with other hematopoietic and non-hematopoietic cell types, suggesting that CD68 may not be a macrophage-specific antigen. CD68 binds to tissue- and organ-specific lectins or selectins, allowing homing of macrophage subsets to particular sites. Rapid recirculation of CD68 from endosomes and lysosomes to the plasma membrane may allow macrophages to crawl over selectin-bearing substrates or other cells.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481