

Anti-ST6GAL1 Antibody (6I659)

Product Details

Ig Type:	Mouse IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	6I659
Purification:	Protein A

Applications

Application:	ELISA
Recommended	ELISA: 1:1000-1:2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human ST6GAL1 / ST6GAL-I protein (TMPY-02093)
Antigen Species:	Human
Synonyms:	ST6 β -galactosamide α -2,6-sialyltransferase 1; ST6 beta-galactosamide α -2,6-sialyltransferase 1; ST6 beta-galactosamide alpha-2,6-sialyltransferase 1

Research Background

Beta-galactoside alpha-2,6-sialyltransferase 1, also known as B-cell antigen CD75, Sialyltransferase 1, CMP-N-acetylneuraminic acid-beta-galactosamide-alpha-2,6-sialyltransferase 1, ST6GAL1 and SIAT1, is a single-pass type II membrane protein that belongs to the glycosyltransferase 29 family. Sialyltransferases are key enzymes in the biosynthesis of sialoglycoconjugates that catalyze the transfer of sialic residue from its activated form to an oligosaccharidic acceptor. ST6GAL1 / SIAT1 is normally found in the Golgi but can be proteolytically processed to a soluble form. It is involved in the generation of the cell-surface carbohydrate determinants and differentiation antigens HB-6, CDw75, and CD76. β -Galactoside α 2,6-sialyltransferases ST6GAL1 and ST6GAL2 are the two unique members of the ST6GAL family described in higher vertebrates. ST6GAL1 / SIAT1 transfers sialic acid from the donor of substrate CMP-sialic acid to galactose containing acceptor substrates.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481