

Anti-HIV-1 (group M, subtype CRF07_BC) gp140 Antibody (8E659)

Product Details

Ig Type:	Mouse IgG1
Conjugation:	Unconjugated
Clone:	8E659
Purification:	Protein A

Applications

Application:	ELISA
Recommended	ELISA: 1:1000-1:2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: HIV-1 gp140 (group M, subtype CRF07_BC) Protein (TMPY-02882)
Antigen Species:	HIV

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481