

Anti-HIV-1 (group M, subtype B, strain 92418) p24 Protein Antibody (2G914)

Product Details

Ig Type:	Mouse IgG1
Conjugation:	Unconjugated
Clone:	2G914
Purification:	Protein A

Applications

Application:	ELISA(Cap)
--------------	------------

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: HIV-1 p24 (group M, subtype B, strain 92418) Protein (TMPY-01552)
Antigen Species:	HIV

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481