

## Anti-HIV-1 (group M, subtype B, strain 92418) p24 Protein Antibody-Biotin (8F825)

### Product Details

Ig Type:	Mouse IgG1
Conjugation:	Biotin
Clone:	8F825
Purification:	Protein A

### Applications

Application:	ELISA
Recommended	ELISA: 0.1-1 µg/ml

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: HIV-1 p24 (group M, subtype B, strain 92418) Protein (TMPY-01552)
Antigen Species:	HIV

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481