

Anti-HIV-1 (group M, subtype B, strain 92418) p24 Protein Antibody-HRP (8F825)

Product Details

Ig Type: Mouse IgG1
Conjugation: HRP
Clone: 8F825
Purification: Protein A

Applications

Application: ELISA
Recommended: ELISA: 0.1-1 µg/ml

Properties

Stability & Storage: Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free. Keep away from direct sunlight.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen: Recombinant Protein: HIV-1 p24 (group M, subtype B, strain 92418) Protein (TMPY-01552)
Antigen Species: HIV

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481