

Anti-HIV-1 (group M, subtype B, strain 92418) p24 Antibody (7N248)

Product Details

Ig Type:	Rabbit IgG
Conjugation:	Unconjugated
Clone:	7N248
Purification:	Protein A

Applications

Application:	ELISA, WB
Recommended	WB: his antibody can be used at 1:500-1:1000 with the appropriate secondary reagents to detect HIV-1 Gag-p24 in WB. Using a DAB detection system, the detection limit for HIV-1 Gag-p24 is approximately 20 ng/lane under non-reducing conditions and 5 ng/lane under reducing conditions.; ELISA: 1:5000-1:10000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: HIV-1 Gag-P24 protein (TMPY-01552)
Antigen Species:	HIV

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481