

Anti-VNN2 Antibody-PE (9J30)

Product Details

Ig Type:	Mouse IgG1
Reactivity:	Human
Conjugation:	PE
Clone:	9J30
Purification:	Protein A

Applications

Verified Activity:	Flow cytometric analysis of anti-human VNN2 on human whole blood monocytes (left panel) and granulocytes (right panel). The fluorescence histograms were derived from events with the forward and side light-scatter characteristics of viable monocytes and granulocytes.
Application:	FCM
Recommended	5 µl/Test, 0.1 mg/ml

Properties

Stability & Storage:	Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight. Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human VNN2 / Vanin-2 protein (TMPY-01614)
Antigen Species:	Human
Synonyms:	VNN2;FOAP-4;vanin 2;GPI-80

Research Background

Vascular non-inflammatory molecule 2, also known as glycosyl-phosphatidyl inositol-anchored protein GPI-8, Vanin-2, Protein FOAP-4 and VNN2, is a cell membrane protein that belongs to the CN hydrolase family and Vanin subfamily. VNN2 is widely expressed with higher expression in spleen and blood. VNN2 is a member of the vanin family of proteins which share extensive sequence similarity with each other, and also with biotinidase. The family includes secreted and membrane-associated proteins, a few of which have been reported to participate in hematopoietic cell trafficking. No biotinidase activity has been demonstrated for any of the vanin proteins, however, they possess pantetheinase activity, which may play a role in oxidative-stress response. VNN2 is an amidohydrolase that hydrolyzes specifically one of the carboamide linkages in D-pantetheine thus recycling pantothenic acid (vitamin B5) and releasing cysteamine. It is involved in the thymus homing of bone marrow cells. VNN2 plays a role in transendothelial migration of neutrophils and may regulate beta-2 integrin-mediated cell adhesion, migration and motility of neutrophil.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481