

Anti-HAI-1/SPINT1 Antibody (7G234)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	7G234
Purification:	Protein A

Applications

Application:	ELISA
Recommended	ELISA: 1:5000-1:10000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human HAI-1 / SPINT1 protein (TMPY-01917)
Antigen Species:	Human
Synonyms:	serine peptidase inhibitor, Kunitz type 1;MANSC2;HAI;SPINT1;UNQ223/PRO256;HAI1
Biology Area:	Serine Proteases and Regulators

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481