

Anti-HAI-1/SPINT1 Antibody (5J371)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	5J371
Purification:	Protein A

Applications

Application:	ELISA(Det)
--------------	------------

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human HAI-1/SPINT1 Protein (TMPY-01917)
Antigen Species:	Human
Synonyms:	HAI1;MANSC2;HAI;serine peptidase inhibitor, Kunitz type 1;SPINT1;UNQ223/PRO256
Biology Area:	Serine Proteases and Regulators

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481