

Anti-Siglec-5 Antibody-FITC (3K993)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	FITC
Clone:	3K993
Purification:	Protein A

Applications

Verified Activity:	Flow cytometric analysis of Human SIGLEC5(CD170) expression on human whole blood granulocytes. Cells were stained with FITC-conjugated anti-Human SIGLEC5(CD170). The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of viable granulocytes.
Application:	FCM
Recommended	5 µl/Test, 0.1 mg/ml

Properties

Stability & Storage:	Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight. Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human SIGLEC5 Protein (TMPY-03549)
Antigen Species:	Human
Synonyms:	OBBP2;SIGLEC-5;CD170;CD33L2;OBBP;SIGLEC5;sialic acid binding Ig like lectin 5;OB-BP2
Biology Area:	ITIM/ITAM Immunoreceptors and Related Molecules

Research Background

SIGLEC5 contains 2 Ig-like C2-type (immunoglobulin-like) domains and 1 Ig-like V-type (immunoglobulin-like) domain. It belongs to the immunoglobulin superfamily and SIGLEC (sialic acid binding Ig-like lectin) family. SIGLEC5 is expressed by monocytic/myeloid lineage cells. It is found at high levels in peripheral blood leukocytes, spleen, bone marrow and at lower levels in lymph node, lung, appendix, placenta, pancreas and thymus. It is also expressed by monocytes and neutrophils but absent from leukemic cell lines representing early stages of myelomonocytic differentiation. SIGLEC5 is a putative adhesion molecule that mediates sialic-acid dependent binding to cells. It binds equally to alpha-2,3-linked and alpha-2,6-linked sialic acid. The sialic acid recognition site may be masked by cis interactions with sialic acids on the same cell surface.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481