

Anti-LYPD3 Antibody-PE (5Y943)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	PE
Clone:	5Y943
Purification:	Protein A

Applications

Verified Activity:	Flow cytometric analysis of Human LYPD3 expression on PC-3 cells.
Application:	FCM
Recommended	10 µl/Test, 0.1 mg/ml

Properties

Stability & Storage:	Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human C4.4A / LYPD3 protein (TMPY-01859)
Antigen Species:	Human
Synonyms:	C4.4A;LY6/PLAUR domain containing 3

Research Background

Ly6 / PLAUR domain-containing protein 3, also known as GPI-anchored metastasis-associated protein C4.4A homolog, Matrigel-induced gene C4 protein, MIG-C4, and LYPD3, is a cell membrane protein that contains two UPAR/Ly6 domains. Human LYPD3 contains two UPAR/Ly6 domains. LYPD3 is expressed in the placenta, skin, and urothelium. It is found in suprabasal keratinocytes of chronic wounds. Weak expression of LYPD3 is found in the esophagus and peripheral blood mononuclear cells. It is found in the majority of primary and metastatic transitional cell carcinomas (TCCs) and as well in breast cancer tissues, but not in adjacent normal tissues. High expression of LYPD3 is found in the tumor component of some noninvasive superficial lesions and invasive and metastatic urothelial cancers. LYPD3 is up-regulated in migrating keratinocytes during epithelisation of incisional skin wounds. LYPD3 supports cell migration. It may be involved in urothelial cell-matrix interactions. It may also be involved in tumor progression

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481