

Anti-Caspase-14 Antibody (1S873)

Product Details

| | |
|---------------|--------------|
| Ig Type: | Mouse IgG2b |
| Reactivity: | Human |
| Conjugation: | Unconjugated |
| Clone: | 1S873 |
| Purification: | Protein A |

Applications

| | |
|--------------|----------------------|
| Application: | ELISA |
| Recommended | ELISA: 1:1000-1:2000 |

Properties

| | |
|----------------------|--|
| Stability & Storage: | Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free. |
| Shipping: | Shipping with blue ice. |

Antigen Details

| | |
|------------------|---|
| Immunogen: | Recombinant Protein: Human Caspase-14 / CASP14 protein (TMPY-01824) |
| Antigen Species: | Human |
| Synonyms: | caspase 14, apoptosis-related cysteine peptidase;ARCI12 |

Research Background

Caspase 14 is a member of the caspase family. Caspases are a kind of cysteine proteinase consisting of a prodomain plus large and small catalytic subunits, that play a central role in cell apoptosis. Caspase 14 possesses an unusually short prodomain and is highly expressed in embryonic tissues but absent from most of the adult tissues except for the skin, which suggests a role in ontogenesis and skin physiology. Unlike the other short prodomain caspases (caspase-3, caspase-6, and caspase-7), Caspase 14 was not processed by multiple death stimuli including activation of members of the tumor necrosis factor receptor family and expression of proapoptotic members of the bcl-2 family. Caspase 14 has been described to be processed and activated by anti-Fas agonist antibody or TNF-related apoptosis inducing ligand in vivo. The expression and processing of this caspase may take part in keratinocyte terminal differentiation, which is essential for the skin barrier.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481