

## Anti-B2M/beta 2-Microglobulin Antibody (4U533)

## Product Details

Ig Type:	Mouse IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	4U533
Purification:	Protein A

## Applications

Application:	ELISA(Cap)
--------------	------------

## Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

## Antigen Details

Immunogen:	Recombinant Protein: Human Beta-2 microglobulin/B2M Protein (TMPY-01686)
Antigen Species:	Human
Synonyms:	beta-2 microglobulin;β-2 microglobulin

## Research Background

B2M, also known as β2-Microglobulin or CDABP0092, is a component of MHC class I molecules found expression in all nucleated cells (excludes red blood cells). The major function of MHC class I molecules is to display fragments of proteins from within the cell to T-cells and cells containing foreign proteins will be attacked. B2M (β2-Microglobulin) is a low molecular weight protein. It was demonstrated that B2M (β2-Microglobulin) was localized in the membranes of nucleated cells and was found to be associated with HL-A antigens. B2M (β2-Microglobulin) is present in free form in various body fluids and as a subunit of histocompatibility antigens on cell surfaces lateral to the α3 chain. Unlike α3, β2 has no transmembrane region. Directly above β2 lies the α1 chain, which itself is lateral to the α2. In the absence of B2M (β2 microglobulin), very limited amounts of MHC class I (classical and non-classical) molecules can be detected on the surface. In the absence of MHC class I, CD8 T cells, a subset of T cells involved in the development of acquired immunity cannot develop. Low levels of B2M (β2 microglobulin) can indicate non-progression of HIV.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481