

## Anti-SMOC1 Antibody (5Z35)

## Product Details

Ig Type:	Mouse IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	5Z35
Purification:	Protein A

## Applications

Application:	ELISA
Recommended	ELISA: 1:1000-1:2000

## Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

## Antigen Details

Immunogen:	Recombinant Protein: Human SMOC-1 Protein (TMPY-01783)
Antigen Species:	Human
Synonyms:	OAS;SPARC related modular calcium binding 1
Biology Area:	Calcium-binding Proteins and Related Molecules

## Research Background

SPARC-related modular calcium-binding protein 1, also known as secreted modular calcium-binding protein 1 and SMOC1, is a member of the SPARC family. SMOC1 is widely expressed in many tissues with a strongest signal in ovary. It contains two EF-hand domains, one Kazal-like domain and two thyroglobulin type-1 domains. Extracellular matrix proteins have been implicated in the regulation of osteoblast differentiation of bone marrow-derived mesenchymal stem cells (BMSCs) through paracrine or autocrine mechanisms. SMOC1 is a regulator of osteoblast differentiation of BMSCs. SMOC1 is highly expressed and secreted in BMSCs stimulated with osteogenic medium (OSM). SMOC1 and SMOC2 are matricellular proteins thought to influence growth factor signaling, migration, proliferation, and angiogenesis. SMOC1 and SMOC2 may mediate intercellular signaling and cell type-specific differentiation during gonad and reproductive tract development.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481