

## Anti-DDOST Antibody (5B994)

## Product Details

|               |              |
|---------------|--------------|
| Ig Type:      | Mouse IgG2a  |
| Reactivity:   | Human        |
| Conjugation:  | Unconjugated |
| Clone:        | 5B994        |
| Purification: | Protein A    |

## Applications

|              |                      |
|--------------|----------------------|
| Application: | ELISA                |
| Recommended  | ELISA: 1:1000-1:2000 |

## Properties

|                      |  |
|----------------------|--|
| Stability & Storage: | Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free. |
| Shipping:            | Shipping with blue ice.  |

## Antigen Details

|                  |  |
|------------------|--|
| Immunogen:       | Recombinant Protein: Human DDOST / OST48 protein (TMPY-02050)  |
| Antigen Species: | Human  |
| Synonyms:        | OST48;CDG1R;WBP1;OKSWcl45;AGER1;dolichyl-diphosphooligosaccharide--protein glycosyltransferase subunit (non-catalytic);OST |

## Research Background

The enzyme oligosaccharyltransferase (dolichyl-diphosphooligosaccharide-protein glycosyltransferase) (DDOST), or 48-kDa subunit (OST48) is one of the catalytic subunits in this complex, exerts a typical type I membrane topology, containing a large luminal domain, a hydrophobic transmembrane domain and a short cytosolic peptide tail. DDOST/OST48 catalyzes the transfer of a high-mannose oligosaccharide (GlcNac2Man9Glc3) from a dolichol-linked oligosaccharide donor (dolichol-P-GlcNac2Man9Glc3) onto the asparagine acceptor site within an Asn-X-Ser/Thr consensus motif in nascent polypeptide chains across the membrane of the endoplasmic reticulum. The mammalian oligosaccharyltransferase (OST) is an oligomeric complex composed of three type I transmembrane proteins of the endoplasmic reticulum: ribophorin I (RI), ribophorin II (RII), and OST48. OST48 is not a glycoprotein and is not recognized by antibodies to either ribophorin. Like ribophorins I and II, OST48 was found to be an integral membrane protein, with the majority of the polypeptide located within the lumen of the endoplasmic reticulum (ER). OST48 does not show significant amino acid sequence homology to either ribophorin I or II. It had been found that only the luminal domain of RI contains ER retention information. The dilysine motif in OST48 functions as an ER localization motif because OST48 in which the two lysine residues are replaced by serine (OST48ss) is no longer retained in the ER and is found instead also at the plasma membrane.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481