

Anti-Ferritin light chain Antibody (4F677)

Product Details

Ig Type:	Mouse IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	4F677
Purification:	Protein A

Applications

Application:	ELISA
Recommended	ELISA: 1:1000-1:2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human FTL / ferritin, light polypeptide protein (TMPY-02582)
Antigen Species:	Human
Synonyms:	LFTD;MGC71996;NBIA3;ferritin, light polypeptide

Research Background

Ferritin, light polypeptide (FTL) is the light subunit of the ferritin protein. Ferritin is the major intracellular iron storage protein in prokaryotes and eukaryotes. It is composed of 24 subunits of the heavy and light ferritin chains. Storage of iron in the tissues occurs in the form of ferritin and hemosiderin. The latter originates from ferritin that has undergone intracellular digestion of its protein shell, leaving the iron core. Ferritin and hemosiderin are components of a continuum. Ferritin has been identified in all types of living organisms: animals, plants, molds, and bacteria. Within the protein shell of ferritin, iron is first oxidized to the ferric state for storage as ferric oxyhydroxide. Thus, ferritin removes excess iron from the cell sap where it could otherwise participate in peroxidation mechanisms.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481