

Anti-Thioredoxin 2/TRX2 Antibody (3V817)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	3V817
Purification:	Protein A

Applications

Application:	ELISA
Recommended	ELISA: 1:5000-1:10000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human Thioredoxin-2 / TXN2 protein (TMPY-02084)
Antigen Species:	Human
Synonyms:	MTRX;RP5-1119A7. 13-005;thioredoxin 2;MT-TRX;TRX2

Research Background

Thioredoxin-2, also known as TXN2, MTRX and TRX2, is a member of the thioredoxin family. Tryparedoxins (TXN) are thioredoxin-related proteins which, as trypanothione:peroxiredoxin oxidoreductases, constitute the trypanothione-dependent antioxidant defense and may also serve as substrates for ribonucleotide reductase in trypanosomatids. Thioredoxin-2 / TXN2 contains one thioredoxin domain. It is widely expressed in adult (at protein level) and fetal tissues. Human Thioredoxin-2 / TXN2 is a small redox protein important in cellular antioxidant defenses, as well as in the regulation of apoptosis. Thioredoxin-2 / TXN2 has an anti-apoptotic function and plays an important role in the regulation of mitochondrial membrane potential. Thioredoxin-2 / TXN2 could be involved in the resistance to anti-tumor agents. It possesses a dithiol-reducing activity. Thioredoxin-2 / TXN2 plays an important role in protecting the mitochondria against oxidative stress and in sensitizing the cells to ROS-induced apoptosis. Mammalian Thioredoxin-2 / TXN2 is a mitochondrial isoform of highly evolutionary conserved thioredoxins. Thioredoxins are small ubiquitous protein-disulfide oxidoreductases implicated in a large variety of biological functions.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481