

Anti-Ezrin Antibody (7L736)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	7L736
Purification:	Protein A

Applications

Verified Activity:	Immunochemical staining of human EZR in human brain with rabbit monoclonal antibody (1:200, formalin-fixed paraffin embedded sections).
Application:	IHC-P
Recommended	IHC-P: 1:100-1:500

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human Ezrin / villin 2 Protein (TMPY-03053)
Antigen Species:	Human
Synonyms:	ezrin;CVIL;CVL;HEL-S-105;VIL2

Research Background

VIL2, also known as villin-2, belongs to the ERM protein family. It is expressed in cerebral cortex, basal ganglia, hippocampus, hypophysis, and optic nerve. VIL2 contains 1 FERM domain and may be involved in connections of major cytoskel

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481