

Anti-Serpin B4 Antibody (7X812)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	7X812
Purification:	Protein A

Applications

Verified Activity:	<ol style="list-style-type: none">1. Immunochemical staining of human SERPINB4 in human skin with rabbit monoclonal antibody (1:200, formalin-fixed paraffin embedded sections).2. Immunochemical staining of human SERPINB4 in human esophageal carcinoma with rabbit monoclonal antibody (1:200, formalin-fixed paraffin embedded sections).
Application:	IHC-P
Recommended	IHC-P: 1:100-1:500

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human SerpinB4 / SCCA2 protein (TMPY-02425)
Antigen Species:	Human
Synonyms:	LEUPIN;SCCA1;SCCA-2;PI11;serpin peptidase inhibitor, clade B (ovalbumin), member 4;Serpin B4;SCCA2

Research Background

Overexpression of SERPINB4 in tumor cells inhibited recombinant GrM-induced as well as NK cell-mediated cell death and this inhibition depended on the reactive center loop of the serpin.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481