

## Anti-IL-20 Antibody (6I697)

### Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	6I697
Purification:	Protein A

### Applications

Application:	ELISA
Recommended	ELISA: 1:5000-1:10000

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Human IL-20 / Interleukin-20 Protein (TMPY-02959)
Antigen Species:	Human
Synonyms:	interleukin 20

### Research Background

IL20/Interleukin-20 belongs to the IL-10 family. It is a cytokine structurally related to interleukin 10. IL20/Interleukin-20 can be detected in skin, trachea, and other tissues. It is produced by activated keratinocytes and monocytes and transmits an intracellular signal through two distinct cell-surface receptor complexes on keratinocytes and other epithelial cells. It has been shown that interleukin-20 transduces its signal through signal transducer and activator of transcription 3 (STAT3) in keratinocytes. It may be involved in epidermal function and psoriasis. It also regulates the proliferation and differentiation of keratinocytes during inflammation, particularly inflammation associated with the skin. Also, IL20/Interleukin-20 also causes cell expansion of multipotential hematopoietic progenitor cells. A specific receptor for this cytokine is found to be expressed in skin and upregulated dramatically in psoriatic skin, suggesting a role for this protein in epidermal function and psoriasis.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481