

Anti-GOLPH2/GOLM1 Antibody-FITC (9P207)

Product Details

| | |
|---------------|-------------|
| Ig Type: | Mouse IgG2a |
| Reactivity: | Human |
| Conjugation: | FITC |
| Clone: | 9P207 |
| Purification: | Protein A |

Applications

| | |
|--------------------|--|
| Verified Activity: | Flow cytometric analysis of Human GOLM1/GP73 expression in HeLa cells. The cells were treated according to manufacturer's manual (BD Pharmingen™ Cat. No. 554714), and then stained with FITC Mouse anti-Human GOLM1/GP73. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of intact cells. |
| Application: | FCM |
| Recommended | 10 µl/Test, 0.1 mg/ml |

Properties

| | |
|----------------------|--|
| Stability & Storage: | Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight. Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal. |
| Shipping: | Shipping with blue ice. |

Antigen Details

| | |
|------------------|--|
| Immunogen: | Recombinant Protein: Human GOLM1 / GP73 protein (TMPY-02193) |
| Antigen Species: | Human |
| Synonyms: | GOLPH2;bA379P1.3;golgi membrane protein 1;C9orf155;HEL46;GP73;PSEC0257 |

Research Background

Golgi membrane protein 1, also known as Golgi membrane protein GP73, Golgi phosphoprotein 2, and GOLM1, is a protein that belongs to the GOLM1 / CASC4 family. GOLM1 is widely expressed. It is highly expressed in the colon, prostate, trachea, and stomach. It is expressed at a lower level in testis, muscle, lymphoid tissues, white blood cells, and spleen. It is predominantly expressed by cells of the epithelial lineage. GOLM1 is expressed at a low level in the normal liver. Expression significantly increases in virus (HBV, HCV) infected liver. Expression of GOLM1 does not increase in liver disease due to non-viral causes (alcohol-induced liver disease, autoimmune hepatitis). Increased expression in hepatocytes appears to be a general feature of advanced liver disease. In liver tissue from patients with adult giant-cell hepatitis (GCH), GOLM1 is strongly expressed in hepatocyte-derived syncytial giant cells. GOLM1 is constitutively expressed by biliary epithelial cells but not by hepatocytes.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481