

Anti-GOLPH2/GOLM1 Antibody (3U78)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	3U78
Purification:	Protein A

Applications

Verified Activity:	<p>1. Immunochemical staining of human GOLM1 in human uterus with rabbit monoclonal antibody (1:200, formalin-fixed paraffin embedded sections). Positive staining was localized to cytoplasm of uterine gland. The left panel: tissue incubated with primary antibody; The right panel: tissue incubated with the mixture of primary antibody and antigen (recombinant protein).</p> <p>2. Immunochemical staining of human GOLM1 in human liver with rabbit monoclonal antibody (1:200, formalin-fixed paraffin embedded sections). Positive staining was localized to cytoplasm of hepatocyte and interlobular bile duct. The left panel: tissue incubated with primary antibody; The right panel: tissue incubated with the mixture of primary antibody and antigen (recombinant protein).</p>
Application:	IHC-P
Recommended	IHC-P: 1:100-1:500

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human GOLM1 / GP73 protein (TMPY-02193)
Antigen Species:	Human
Synonyms:	golgi membrane protein 1;GP73;C9orf155;GOLPH2;bA379P1.3;PSEC0257;HEL46

Research Background

Golgi membrane protein 1, also known as Golgi membrane protein GP73, Golgi phosphoprotein 2, and GOLM1, is a protein that belongs to the GOLM1 / CASC4 family. GOLM1 is widely expressed. It is highly expressed in the colon, prostate, trachea, and stomach. It is expressed at a lower level in testis, muscle, lymphoid tissues, white blood cells, and spleen. It is predominantly expressed by cells of the epithelial lineage. GOLM1 is expressed at a low level in the normal liver. Expression significantly increases in virus (HBV, HCV) infected liver. Expression of GOLM1 does not increase in liver disease due to non-viral causes (alcohol-induced liver disease, autoimmune hepatitis). Increased expression in hepatocytes appears to be a general feature of advanced liver disease. In liver tissue from patients with adult giant-cell hepatitis (GCH), GOLM1 is strongly expressed in hepatocyte-derived syncytial giant cells. GOLM1 is constitutively expressed by biliary epithelial cells but not by hepatocytes.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481