

Anti-VSTM1 Antibody (6K35)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	6K35
Purification:	Protein A

Applications

Application:	ELISA
Recommended	ELISA: 1:5000-1:10000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human VSTM1 protein (TMPY-02555)
Antigen Species:	Human
Synonyms:	SIRL1;V-set and transmembrane domain containing 1;SIRL-1;VSTM1;UNQ3033;MGC119161;MGC119160

Research Background

V-set and transmembrane domain containing 1 (VSTM1) is a protein containing the V-set domains. V-set domains are immunoglobulin-like domains resembling the antibody variable region. V-set domains are found in many kinds of protein families, including immunoglobulin light and heavy chains, several T-cells such as CD2, CD4, CD80, and CD86, myelin membrane adhesion molecules, junction adhesion molecules (JAM), tyrosine-protein kinase receptors, and the programmed cell death protein1.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481