

Anti-Thyroid peroxidase Antibody (7P893)

Product Details

Ig Type:	Mouse IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	7P893
Purification:	Protein A

Applications

Verified Activity:	Immunochemical staining of human TPO in human thyroid gland with mouse monoclonal antibody (1:60, formalin-fixed paraffin embedded sections).
Application:	ELISA(Cap),IHC-P
Recommended	IHC-P: 1:50-1:200; ELISA(Cap): 1:250-1:2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human Thyroid peroxidase / TPO Protein (TMPY-02779)
Antigen Species:	Human
Synonyms:	MSA;TPX;thyroid peroxidase;TDH2A

Research Background

Thyroid peroxidase is a membrane-bound glycoprotein which belongs to the peroxidase family, XPO subfamily. It contains 1 EGF-like domain and 1 Sushi (CCP/SCR) domain. Thyroid Peroxidase represents one of the main autoantigenic targets in autoimmune thyroid disease of humans. It used to be taken as the formerly so-called `microsomal antigen` several years ago. As an integral membrane glycoprotein it is restricted to the apical plasma membrane of the follicular epithelial cells and comprises two identical subunits of approx 100 kDa molecular weight. Thyroid peroxidase is an enzyme expressed abundantly in the thyroid that liberates iodine for addition onto tyrosine residues on thyroglobulin for the production of thyroxine or triiodothyronine, thyroid hormones. Thyroid peroxidase plays a key role in the thyroid hormone biosynthesis by catalysing both the iodination of tyrosyl residues and the coupling of iodotyrosyl residues in thyroglobulin to form precursors of the thyroid hormones T4 and T3.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481